

231/3 – BIOLOGY FORM THREE

CONFIDENTIAL

Moss plant labeled K in a watch glass / peri-dish. Showing (Rhizoids, leaves, seta, capsule.

Hand lens

REQUIREMENTS

Each candidate will require the following:

- 5 test tubes
- 5ml of solution L₁ starch solution 0.5%
- 5ml of solution L₂ diastase solution 10%
- 2ml of solution L₃ Boiled diastase solution 10%

NB: boiled for 10 minutes

- Means of labeling 3 labels
- Water bath
- Thermometer
- Iodine solution
- 3 Droppers
- A White tile
- Means of timing